



prepFAST M5

Syringe Driven Autodilution




prepFAST M5 Inline Autodilution

prepFAST M5 Automation


The prepFAST M5 has an additional syringe for sample loading allowing for precise and accurate loading of micro volume and/or viscous samples. The 5th syringe also provide additional features such as premixing samples and enhanced washout for lower blanks. Combined with the autocalibration, inline sample dilution and automatic overrange dilution the capabilities of prepFAST, the prepFAST M5 is an advanced laboratory automation system.

FAST Uptake & Washout

Normal Analysis



FAST Analysis



Six Steps in a Standard Analysis

1. Autosampler Movement
2. Uptake
3. Stabilization
4. **Measurement**
5. Wash
6. Overhead

prepFAST has the quickest signal stabilization and most complete rinse-out of any sample introduction system, reducing uptake and wash times while increasing sample throughput.

Benefits

- Fully Automated
 - Autocalibrate
 - Autodilute
 - QC Autodilute
- Rapid uptake, stabilization and washout
- Inline dilution
- Integrated with ICP/ICPMS Software



S500V2 Syringe Pump

The precise ($\pm 0.05\%$) and accurate ($\pm 0.2\%$) S500V2 syringe pump provides smooth and balanced delivery of solutions over a wide range of flow rates (1 to 40,000 $\mu\text{L}/\text{min}$) to ensure rapid reliable inline dilutions.

Five Syringes:

1. Rinse
2. Carrier
3. Diluent
4. Internal Standard
5. Sample Uptake



Syringe Loading

- Premix Sample Before Loading
- Precise Loading of Samples
 - Load Micro Volumes of Sample
 - Operate with Microtiter Plates
 - Load Viscous Samples
- Superior Washout
- Lower Blanks

Syringe Dilution

- Clean
- Chemically Resistant
 - Organics
 - Strong Acids
- Low Maintenance
- Accurate and Precise
- Long-Term Stability

Fluoropolymer P-Series Valves


- PFA rotors for long life and low maintenance
- Chemically resistant compression ring
 - Uniform pressure
 - Keyed for easy assembly
 - Prevents damage from overtightening
- Internal components PTFE coated to prevent corrosion
- All fluoropolymer flow path




Automated Inline Autodilution

Inline Dilutions


1. Syringe load sample into loop




2. Syringes add internal standard, diluent and sample into second loop



3. Diluted sample is injected and sample loop is cleaned




4. Sample and dilution loops are cleaned



*Patent Pending


Intelligent Autodilution

Automated QC Dilutions



Flowchart illustrates the workflow logic for intelligent autodilution of samples that fail QC.

Enhanced Washout



Washout greater than 20,000x is achieved in 10 seconds. This is more than sufficient to eliminate any carry-over from the high concentration QC over range sample onto the intelligently diluted re-run.

Autocalibration

Flexible autocalibration strategies are easily automated and employed with the prepFAST M5.

- Simplifies complex calibration routines such as USP bracketing and multiple stock calibrations
- Only prepares calibration standards as needed – improves stability and reduces lab waste
- On-line dilutions remove risk of contamination from laboratory sources

Conventional Calibration vs prepFAST Autocalibration

Conventional Calibration (seven points)

Offline Prep: 1 Blank + 7 Standards

	STD Position	Concentration
	10	0
	1	2.5
	2	5
	3	12.5
	4	50
	5	100
	6	200
	7	500

prepFAST Autocalibration (seven points)

Inline Prep: 1 Blank + 1 Standard

	STD Position	Inline Dilution Factor	Concentration
	10	200x	0
	2	200x	2.5
	2	100x	5
	2	40x	12.5
	2	10x	50
	2	5x	100
	2	2.5x	200
	2	1x	500

prepFAST Calibration: Blank can be analyzed diluted or undiluted.

Multiple Curves from a Single Stock Standard

USP <233> Conventional Calibration (4 Dosage Levels)

Offline Prep: 1 Blank + 8 Standards

Daily Dose	STD Position	STD	As Concentration (µg/L)
	10		0
20g/day	1	0.5J	0.075
	2	2J	0.3
10g/day	3	0.5J	0.15
	4	2J	0.6
5g/day	5	0.5J	0.3
	6	2J	1.2
1g/day	7	0.5J	1.5
	8	2J	6

USP <233> prepFAST Autocalibration (4 Dosage Levels)

Inline Prep: 1 Blank + 1 Standard


Daily Dose	STD Position	Inline Dilution Factor	As Concentration (µg/L)
	10	1x	0
20g/day	1	80x	0.075
	1	20x	0.3
10g/day	1	1x	0
	1	40x	0.15
	1	10x	0.6
5g/day	1	1x	0
	1	20x	0.3
	1	5x	1.2
1g/day	1	1x	0
	1	4x	1.5
	1	1x	6

prepFAST Calibration: Blank can be analyzed diluted or undiluted.

Calibration Using Multiple Standards


Conventional Calibration (seven points)

Offline Prep: 1 Blank + 7 Standards




prepFAST Autocalibration (seven points)

Inline Prep: 1 Blank + 1 Standard



prepFAST Calibration: Blank can be analyzed diluted or undiluted.


Autocalibration ²³⁸U



Predefined dilution factors for a single multi-element standard are used to build calibration curves. Accuracy of dilution results in high linearity (R = 1.0000).


Environmental

Trace Metals in Water CRM




Precise and accurate results for concentrations of metals in a Trace Metals in Water Certified Reference Material (CRM, NIST 1643e) are obtained with prepFAST autocalibration and autodilution.

Accurate Spike Recovery




Omaha tap water is used as the EPA 200.8 defined, Laboratory Fortified Matrix (LFM) and Quality Control Standard (QCS). Recovery of 95-105% greatly exceeds EPA requirements of 70-130%.

Sample Stability



Fire assay digests for precious metals in Aqua Regia require dilution before analysis. prepFAST inline 20x autodilution eliminates the manual dilution step.

Accuracy





High TDS ore samples (H_2SO_4 , HNO_3 , HF and H_3BO_3 digests) are placed directly in the prepFAST M5 and autodiluted (10x) at time of analysis on ICP. High recoveries are obtained while improving sample throughput and minimizing operator handling of dangerous chemicals.


Clinical (Urine)

- Automated calibration routine with prepFAST M5 saves significant operator time
- Excellent accuracy and reproducibility demonstrated for urine samples
- Samples diluted on-line by 10x, reduces sample handling

Autocalibration with 100µL Sample Loop



Recovery Using 100µL Urine Samples




Clinical (Blood)

- Automated calibration and sample handling demonstrate excellent accuracy on sample spikes and reference materials
- Blood samples diluted 50x online; reduced sample handling also improves operator safety
- Automated matrix-matched calibration performed on-line with prepFAST M5

Spike Recovery from Undiluted Blood

Spike Blood Sample Recovery Chart



UTAK Blood CRM Accuracy and Precision

UTAK Control Level 2

Sample Id	Cr 52 (ppb)	Mn 55 (ppb)	As 75 (ppb)	Se 78 (ppb)	Se 82 (ppb)	Cd 111 (ppb)	Hg 202 (ppb)	Tl 205 (ppb)	Pb 208 (ppb)
Expected Value (ppb)	19.300	30.400	27.600	207.000	207.000	5.300	15.800	11.500	402.000
% Recovery	94%	95%	103%	100%	107%	96%	102%	95%	99%
% RSD	0.8%	1.4%	4.3%	1.8%	2.1%	2.3%	2.2%	1.0%	1.6%

UTAK Control Level 3

Sample Id	Cr 52 (ppb)	Mn 55 (ppb)	As 75 (ppb)	Se 78 (ppb)	Se 82 (ppb)	Cd 111 (ppb)	Hg 202 (ppb)	Tl 205 (ppb)	Pb 208 (ppb)
Expected Value	69.000	90.000	84.500	604.000	604.000	14.300	51.700	45.200	597.000
% Recovery	95%	98%	101%	104%	106%	95%	98%	94%	99%
% RSD	0.3%	0.9%	1.7%	0.6%	0.3%	1.6%	2.6%	2.4%	3.6%

Pharmaceutical

- prepFAST M5 performs required USP bracketing calibrations standards for each sample dosage
- Provides automated sample dilutions

Autocalibration

- **Calibration** Blank, 0.5J and 1.5J

Automated Inline Sample Dilution and J Spiking


- **Suitability** Drift determined at 1.5J
 - NMT 20% after sample analysis
- **Repeatability** determined on six 1J spiked sample replicates
 - NMT 20%
- **Ruggedness** repeatability determined on 2 separate days, operators, instruments
 - NMT 25%
- **Accuracy** determined on six 1J spiked samples
 - NMT 70%-150%
- **Detectability** determined on 0.8J and 1J spiked samples
 - 1J spiked sample must be greater than 0.8J spiked sample



Provides unique calibration protocols for USP elemental impurities analysis.




ICH "J" Calibration Coefficients




Multiple calibration curves for the USP elements have good linearity for drugs with a large range of daily doses. (ICH = International Conference on Harmonisation)



Repeatability Determined on Six 1J Spiked Sample Replicates




Ruggedness Repeatability: 2 Separate Days, Operators, Instruments



Notes

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.



prepFAST M5 Applications

- Environmental
- Pharmaceutical
- Clinical
- Biological
- Geochemical

© Elemental Scientific | 7277 World Communications Drive | Omaha, NE 68122
Tel: 402-991-7800 | sales@icpms.com | www.icpms.com