



Elemental Scientific

ICP | ICPMS | AA

Automated Introduction for PPQ Metals Determination

prepFAST PPQ for Agilent ICPMS

PPQ Level Sample Introduction System for High Purity Chemicals



ICP | ICPMS | AA

prepFAST PPQ

Automated Sample Preparation and Introduction System for PPQ Metals Determination in High Purity Chemicals

The prepFAST PPQ is the most advanced way ultrapure semiconductor grade chemicals are analyzed with ICPMS detection. The prepFAST PPQ utilizes syringe-driven flows of UPW, semiconductor grade acids, and standard solution to automate both sample dilutions and standard curve generation. It eliminates manual handling of samples to deliver sub-ppt detection limit capabilities in direct analysis mode and ppq detection limit capabilities in concentration mode.

Two High Purity Modes

Direct Analysis Mode

- Capability to analyze all semiconductor grade chemicals
- PPT/PPQ detection limits for all semiconductor elements
- General purpose for any sample matrix
- Automated MSA calibration
- Automated inline dilution

Concentration Mode

- PPQ detection limits for analysis of UPW, 30% H₂O₂ and IPA
- Removes difficult matrices such as IPA while keeping metals
- Reduced impact of ICPMS interferences on results
- Improved sensitivity
- Automatic MSA calibration

prepFAST PPQ Autocalibration

The prepFAST PPQ automatically prepares calibration curves for over 40 elements controlled in semiconductor manufacturing processes. Calibrations are generated by automatically diluting an enclosed multielement stock standard. Automation with the high-purity prepFAST PPQ achieves low to sub-ppt calibrations in direct analysis mode and ppq calibrations in concentration mode.



High Purity Automation with PPT/PPQ results

Automation


- Automatic external and MSA calibrations
- Automated sample sensing
 - Accounts for viscosity and automatically adjusts timing
 - Detects and injects the sample and triggers the ICPMS
- Automated syringe-driven sample introduction
 - Sample loading
 - Sample preparation
 - Inline dilution
 - Acid addition (direct mode only)

Ultraclean

- Ultraclean sample preparation
- Integrated ultraclean sample environment
 - Includes ultraclean air shower
 - Options include:
 - ULPA (Ultrapure air) filter
 - Sample racks for PFA containers (<1 mL to 500 mL)
- Continuously-flowing high purity UPW rinse (user-supplied UPW)
- UPW polishing columns for low background

prepFAST PPQ System						
System	Integrated Mobile Autosampler & Enclosure	Ultraclean Air Shower	Integrated FAST valve modules	PFA Nebulizer with Integrated Capillary	PFA Sample Probes	Syringe Pump S500V2
prepFAST PPQ	✓	✓	✓	✓	✓	✓

Pure Automation



prepFAST PPQ Inline Dilution of Semiconductor-grade Chemicals




Diagram showing 10x inline dilution of concentrated H₂SO₄ with prepFAST PPQ.

The prepFAST PPQ allows dilution by volume or weight for IPA and H₂O₂ in concentration mode, and all semiconductor-grade chemicals in direct analysis mode. Metals are quantified using automated inline MSA or external calibration. Automated direct analysis of concentrated chemicals eliminates sample contamination caused by manual dilution into a secondary container.

Examples of Semiconductor Chemicals Analyzed at the ppt Level with prepFAST PPQ*						
Acids	98% H ₂ SO ₄	89% H ₃ PO ₄	70% HNO ₃	49% HF	35% HCl	30% H ₂ O ₂
Bases	22% NH ₄ OH	2.38% TMAH	25%TMAH	KOH		
Organics	IPA	PGMEA/PGME	Photoresist	NMP	Butyl Acetate	Cyclohexanone

*This table contains only a partial list of chemicals which can be analyzed using prepFAST PPQ in direct analysis mode.


Concentration Mode

Example of Detection Limits in Non-cleanroom Environment in Concentration Mode			
Element	DL (ppq)	Element	DL (ppq)
⁷ Li	0.02	⁸⁵ Rb	0.01
¹¹ B	20	⁸⁸ Sr	0.06
²³ Na	0.6	⁹⁰ Zr	0.1
²⁴ Mg	0.2	⁹³ Nb	0.2
²⁷ Al	0.6	⁹⁸ Mo	0.3
²⁸ Si	2700	¹⁰⁷ Ag	0.05
³⁹ K	0.4	¹¹⁴ Cd	0.06
⁴⁰ Ca	0.3	¹¹⁵ In	0.01
⁴⁸ Ti	0.2	¹¹⁸ Sn	0.5
⁵¹ V	0.7	¹²¹ Sb	0.4
⁵² Cr	0.4	¹³³ Cs	0.007
⁵⁵ Mn	0.07	¹³⁸ Ba	0.04
⁵⁶ Fe	0.3	¹⁸⁰ Hf	0.2
⁵⁸ Ni	0.07	¹⁸¹ Ta	0.3
⁵⁹ Co	0.01	¹⁸⁴ W	0.4
⁶³ Cu	0.07	¹⁹⁵ Pt	0.1
⁶⁴ Zn	0.4	²⁰⁵ Tl	0.03
⁶⁹ Ga	0.004	²⁰⁸ Pb	0.05
⁷⁵ As	0.2	²⁰⁹ Bi	0.05

prepFAST PPQ Autocalibration


Autocalibrations for elements controlled in semiconductor manufacturing processes are generated by automatically diluting an enclosed multielement stock standard. Automation with the high-purity prepFAST PPQ achieves ppt to sub-ppt calibration in direct analysis mode and ppq calibrations in concentration mode.

UPW Autocalibration in Concentration Mode




Concentration Mode


UPW Autocalibration in Concentration Mode




UPW Autocalibration in Concentration Mode



UPW Autocalibration for ^{39}K



UPW Autocalibration for ^{58}Ni





Direct Analysis Mode

Example of Detection Limits in Non-cleanroom Environment
in Direct Analysis Mode

Element	DL (ppt)	Element	DL (ppt)
⁷ Li	0.02	⁷² Ge	0.04
⁹ Be	0.007	⁷⁵ As	0.2
¹¹ B	0.9	⁸⁵ Rb	0.008
²³ Na	0.07	⁸⁸ Sr	0.008
²⁴ Mg	0.01	⁹⁰ Zr	0.004
²⁷ Al	0.03	⁹³ Nb	0.002
³⁹ K	0.06	⁹⁵ Mo	0.5
⁴⁰ Ca	0.4	¹¹¹ Cd	0.07
⁴⁸ Ti	0.02	¹¹⁵ In	0.004
⁵¹ V	0.1	¹¹⁸ Sn	0.05
⁵² Cr	0.1	¹²¹ Sb	0.05
⁵⁵ Mn	0.009	¹³⁷ Ba	0.04
⁵⁶ Fe	0.04	¹⁷⁸ Hf	0.003
⁵⁸ Ni	0.01	¹⁸¹ Ta	0.01
⁵⁹ Co	0.007	¹⁸² W	0.01
⁶⁰ Ni	0.01	²⁰⁵ Tl	0.002
⁶³ Cu	0.03	²⁰⁸ Pb	0.005
⁶⁴ Zn	0.04	²³² Th	0.000
⁷¹ Ga	0.002	²³⁸ U	0.003

Autocalibration in Direct Analysis Mode





Elemental Scientific

www.icpms.com

© Elemental Scientific | +1 402-991-7800 | sales@icpms.com | 7277 World Communications Drive | Omaha, NE 68122

19207-4

