



SAMPLE TRAX^{S5}

Barcoded Sample Tracking and
ICPMS Analysis of High Purity
Semiconductor Grade Chemicals



Elemental Scientific

ICP | ICPMS | AA

SAMPLETRAX



SAMPLETRAX S5 Analytical Station

SAMPLETRAX is an advanced, automated sample identification system that uses barcodes to track samples from time of collection through reception to final analysis and data reporting. Direct analysis of semiconductor grade chemicals by ICPMS at less than 1 ppt are achieved with SAMPLETRAX S5.

How it Works

SAMPLETRAX S5 Analytical station automatically:

1. Scans bottles
2. Groups samples by chemical type
3. Analyzes grouped chemicals in a user-defined order
4. Performs wash method specific to each chemical type after each group is analyzed
5. Generates and reports data

Sample Identification

- Barcode scanning accesses information including:
 - Sample type
 - Sample Information (Line / Sampling point / Name / etc.)
 - Method of standardization and analysis

Ultra Pure

- Ultra-clean
- <1 ppt semiconductor metals
- Automated matrix matched MSA, addition or external calibration
- Analytical stations for ICPMS instruments

Laboratory Automation

- Bottle tracking
- Bottle history (cleaning, sample, chemical, analysis, concentration)
- Chemical grouping
- Chemical specific rinse function per chemical
- Customized network, bottle cycle
- Data management



1 mL cup and 250 mL bottle with 2D bottom, 2D cap and 1D side barcode

Fully Automated Sample Identification and Tracking for Ultra-Pure Chemicals




An integrated barcode reader scans the bottom of a PFA bottle to identify sample information before analysis.

Synchronizing the patented barcode reading arm with the sample probe ensures each sample is positively identified at the time of analysis. The enclosed ultra-pure autosampler guarantees sample integrity is not compromised.

Autocalibration of ^{39}K from a Single Stock Solution

SAMPLETRAX S5 analytical station automatically:

1. Scans bottles
2. Groups samples by chemical type
3. Analyzes grouped chemicals in a user-defined order
4. Performs wash method specific to each chemical type after each group is analyzed
5. Generates and reports data



* All calibration strategies automated (MSA, Addition, External)

Barcoded Bottles

A variety of precleaned barcoded bottles, vials and caps are available.



1D and 2D Barcodes Permanently Marked Into Inert PFA

Benefits

- Chemically inert barcodes
- Non-contaminating markings into acid-resistant PFA
- Reusable
- 2D barcoded bottles are compatible with SAMPLETRAX scanning automation systems
- Track bottle position and sample identity
- Precleaned

Types

- Bottles
 - 2D bottom barcoded
 - 1D side barcoded
- Caps
 - Available with 2D barcodes
- Any type of bottle can be custom labeled

TRAX Scan Stations

Scan stations are used to associate sample and analytical information with the bottle's barcode.

Select Sample Type

<input type="checkbox"/>	H2O2
<input checked="" type="checkbox"/>	NH4OH
<input type="checkbox"/>	DSP
<input type="checkbox"/>	LAL
<input type="checkbox"/>	DHF
<input type="checkbox"/>	SC1
<input type="checkbox"/>	IPA
<input type="checkbox"/>	H2SO4
<input type="checkbox"/>	HF

1.) Select sample type



2.) Scan bottle



3.) Add user input

- Location
- Operator
- Bottle Status
- Chemical
- Comments/notes
- Customizable inputs

Multiple scan stations can be used to track bottle usage both in the fab and the laboratory.

TRAX S5 Analytical

Status Indicator Light

Alerts user to system problems
(leaks, empty reagents, etc...)

Integrated Computer Monitor with Swivel Arm

Convenient, protected location
for instrument computer


Storage Drawers

High Purity PFA Reservoirs


With integrated level sensors

Wheels

Allows for flexible positioning



Station




Automated Grouping and Wash


TRAX S5 Analytical Station simplifies and performs the most demanding sample analysis in 3 easy steps.

1. Discover and group samples
2. Create a sequence, a) MSA, b) Addition Cal, c) External Cal.
3. Run samples, acquire data and calculate results

1) Discover Samples



TRAX S5 discovers sample location and analytical information.



Group Samples

Based on discovered sample information, the sequence groups samples by type and sorts by analytical order. Chemical specific washes are performed after each group.


Sample Group Analysis Order	
UPW - Group 1	H ₂ O ₂ - Group 2
2	1
4	3
6	5

MSA Sequence

2a) Select Desired Sequence


External	MSA	Addition	Sequence
----------	-----	----------	----------

Build MSA curves according to protocol for each sample and perform chemical specific washes after each sample type.




Addition Cal Sequence

2b) Select Desired Sequence



Build Addition calibration curves according to protocol for each sample and perform chemical specific washes after each sample type.




External Cal Sequence

2c) Select Desired Sequence Type

External	MSA	Addition	Sequence
----------	-----	----------	----------

Build External calibration curves according to protocol for each sample and perform chemical specific washes after each sample type.



Run Samples


3) Run samples, acquire data and calculate results for reporting to sample

TRAX S5 automatically prepares desired calibration curve and samples for analysis.

As each chemical is prepared and injected, TRAX S5 triggers the 7900 to acquire data and imports raw intensities.

Raw intensities are:

- Associated with TRAX ID
- Used to calculate concentrations
- Stored in database




TRAX database




● Intensity ○ Concentration

Cool NH3...	24Mg-(Cool NH3...	27Al-(Cool NH3...	39K-(Cool NH3...	40Ca-(Cool NH3...	51V-(Hot NH3 1.0...	52Cr-(Cool NH3...	55Mn-(Cool NH3...	56Fe-(Cool NH3...
803	164	155	273	1,182	1,320	43	89	226
863	162	174	292	1,592	1,335	43	94	260
2,813	546	400	681	1,916	1,658	179	283	418
5,258	929	670	1,137	2,540	2,041	351	518	592
11,894	1,993	1,423	2,408	4,429	3,168	807	1,167	1,116
22,860	3,739	2,734	4,435	7,434	4,926	1,611	2,139	1,941
917	103	318	304	1,191	1,315	25	34	298
1,982	272	525	574	12,956	1,354	26	50	380
1,851	165	340	455	4,851	1,264	26	62	348
4,427	219	598	2,874	8,602	1,340	42	39	354
2,105	170	464	687	5,568	1,330	28	38	379
2,214	214	617	460	2,999	1,316	33	38	417
1,098	103	478	340	1,906	1,256	25	33	336



14



Bottle Manager: Tracking

Scan stations update the database in real time allowing the bottle manager to instantaneously provide a snapshot of every bottle's status.

Bottle Tracking

Select "ALL" to view the status of the entire bottle inventory as a color-coded grid.

Quantity of bottles is indicated with (##) in each of the seven color-coded stages.

Bottle History

Click on the TRAXID to view its complete history.

- Status
- Sample type
- Purpose
- Concentrations



12 mL vial



250 mL bottle



1 mL vial

The screenshot shows the 'SAMPLETRAX DataManager' interface. It has a top navigation bar with 'SampleType Editor', 'Receiving', and 'In Analysis'. Below this is a 'Bottle Status' section with a list of color-coded boxes representing different stages: 'ALL (283)' (blue), 'Available (48)' (light blue), 'Carry Out (51)' (dark blue), 'Receiving (61)' (yellow), 'In Analysis (0)' (red), 'Completed (6)' (purple), 'Cleaning (117)' (green), and 'Incomplete (0)' (orange). To the right is a 'Bottle List' section showing a vertical list of bottles with their TRAXIDs: A0280, A0281, MSATEST83, TRAX000027, TRAX000068, TRAX000224, TRAX000135, TRAX000076, TRAX000397, TRAX000305, TRAX000278, TRAX000364, TRAX000349, TRAX000309, and TRAX000328. At the bottom is a 'Page Control' section with navigation buttons and 'Page 1 of 2'.

Clicking a TRAXID opens the bottle history window which shows the complete history of a bottle in the TRAX system including analytical data.

ing and History

Cleaning Available Carry Out DataViewer BottleViewer

A0092	MSATEST32	MSATEST81	TRAX000090	TRAX000006	TRAX000157	TRAX000073	TRAX000162	TRAX000078
AY3AA00029	000002	PFA25B0060	500PF00001	TRAX000057	TRAX000002	TRAX000001	000001	AY3AA00024
PFA25B0001	TRAX000021	TRAX000038	TRAX000026	TRAX000039	TRAX000052	TRAX000019	TRAX000044	TRAX000040
TRAX000029	TRAX000032	TRAX000045	TRAX000053	TRAX000042	TRAX000059	TRAX000050	TRAX000049	TRAX000031
TRAX000054	TRAX000014	TRAX000016	TRAX000051	TRAX000056	TRAX000036	TRAX000023	TRAX000046	TRAX000201
TRAX000196	TRAX000215	TRAX000219	TRAX000218	TRAX000227	TRAX000216	TRAX000197	TRAX000208	TRAX000217
TRAX000091	TRAX000109	TRAX000121	TRAX000115	TRAX000134	TRAX000149	TRAX000123	TRAX000094	TRAX000104
TRAX000074	TRAX000060	TRAX000066	TRAX000041	TRAX000064	TRAX000070	TRAX000063	TRAX000394	TRAX000400
TRAX000399	TRAX000382	TRAX000386	TRAX000376	TRAX000303	TRAX000315	TRAX000321	TRAX000287	TRAX000299
TRAX000294	TRAX000264	TRAX000258	TRAX000253	TRAX000271	TRAX000295	TRAX000301	TRAX000313	TRAX000279
TRAX000274	TRAX000310	TRAX000316	PFA25B0006	TRAX000354	TRAX000366	TRAX000384	TRAX000365	TRAX000347
TRAX000383	TRAX000381	TRAX000363	TRAX000345	TRAX000351	TRAX000380	TRAX000374	TRAX000362	TRAX000338
TRAX000355	TRAX000361	TRAX000367	TRAX000373	TRAX000379	TRAX000385	TRAX000391	TRAX000392	TRAX000297
TRAX000311	TRAX000312	TRAX000270	TRAX000288	TRAX000293	TRAX000292	TRAX000285	TRAX000277	TRAX000329
TRAX000327	TRAX000325	TRAX000273	TRAX000290	TRAX000228	TRAX000177	TRAX000182	TRAX000188	TRAX000213

Bottle History for MSATEST32

TRAX ID	Status	Sample Type	User	Datetime	Protocol	Purpose	Line	Sampling Point	Comments
MSATEST32	Completed	Validation	ESI\Agilent\CPMS	2019-01-29 11:37:44	Semiconductor				
MSATEST32	In Analysis	Validation	ESI\Agilent\CPMS	2019-01-29 11:03:01	Semiconductor				
MSATEST32	Receiving	Validation	ESI\Agilent\CPMS	2019-01-29 10:59:35	Semiconductor				
MSATEST32	Carry Out		ESI\Agilent\CPMS	2019-01-29 10:56:43					
MSATEST32	Available		ESI\Agilent\CPMS	2019-01-29 10:53:07					
MSATEST32	Cleaning		ESI\Agilent\CPMS	2019-01-29 10:46:37					
MSATEST32	Completed	Validation	ESI\Agilent\CPMS	2019-01-29 09:42:04	Semiconductor				
MSATEST32	In Analysis	Validation	ESI\Agilent\CPMS	2019-01-29 09:07:18	Semiconductor				
MSATEST32	Completed	Validation	ESI\Agilent\CPMS	2019-01-28 18:25:44	Semiconductor				
MSATEST32	In Analysis	Validation	ESI\Agilent\CPMS	2019-01-28 17:50:13	Semiconductor				
MSATEST32	In Analysis	Validation	ESI\Agilent\CPMS	2019-01-28 17:47:11	Semiconductor				
MSATEST32	Receiving	Validation	ESI\Agilent\CPMS	2019-01-28 17:37:20	Semiconductor				


TraxId	SampleType	Device	DateTime	Standard	Na	Mg	N	K	Ti	V	Cr	Fe	Ni	Co
MSATEST32	Validation	2	1/28/2019 5:56:42 PM	std-0	0.184	0.001	0.002	0.004	0.002	0.001	0.056	0.002	0	0
MSATEST32	Validation	2	1/29/2019 9:13:03 AM	std-0	0.139	0.108	0.124	0.042	0.108	0.016	0.212	0.113	0.041	0.111

Close

Flexible, Customizable SAM

Automate the whole lab by incorporating

- Network multiple SAMPLETRAX systems
- Network single or multiple labs
- Network SAMPLETRAX with laboratory information management systems (LIMS)
- Standardize sample preparation, calibration and calculation throughout network



MPLETRAX Networking

multiple scanning and analytical stations



Multiple
Analytical Stations



Full System at a Glance

Features and Benefits

Sample Identification

- Barcode scanning accesses information including:
 - Sample type
 - Sample Information (Line, sampling point, name, etc.)
 - Method of standardization and analysis

Ultra Pure

- Ultra-clean
- <1 ppt semiconductor metals
- Automated matrix matched MSA or external calibration
- Analytical stations for ICPMS instruments

Laboratory Automation

- Bottle tracking
- Bottle history (cleaning, sample, chemical, analysis, concentration)
- Chemical grouping
- Chemical specific rinse function per chemical
- Customized network, bottle cycle

How it Works

SAMPLETRAX S5 Analytical station automatically:

1. Scans bottles
2. Groups samples by chemical type
3. Analyzes grouped chemicals in a user-defined order
4. Performs wash method specific to each chemical type after each group is analyzed
5. Generates and reports data

**Examples of Semiconductor Chemicals
Analyzed at the ppt Level with SAMPLETRAX S5**

Acids	Bases	Organics	Chemical Mixes
98% H ₂ SO ₄	22% NH ₄ OH	IPA	SC-1
89% H ₃ PO ₄	2.38% TMAH	PGMEA/PGME	SC-2
70% HNO ₃	25% TMAH	Photoresist	BOE
49% HF	KOH	NMP	DSP
35% HCl		Butyl Acetate	
30% H ₂ O ₂		Cyclohexanone	
SPM			
FPM			
DHF			
Etchant			
Others			

All semiconductor pure chemicals can be analyzed using SAMPLETRAX S5. This table contains only a partial list of common chemicals.

Notes

Notes



Elemental Scientific

ICP | ICPMS | AA

Contact us by phone at 402.991.7800 or by e-mail at sales@icpms.com. Our scientists and engineers are available to answer your questions related to elemental analysis. We are pleased to provide our customers complimentary analytical advice from our on-staff chemists.



© Elemental Scientific | 7277 World Communications Drive | Omaha, NE 68122
Tel: 402-991-7800 | sales@icpms.com | www.icpms.com