

CleanAcids Features

PTFE Elbow

- Connects PFA Source & Product bottles for enclosed environment

Semi-Permeable Membrane

- Releases pressure without venting acid vapor

CleanClever Well

- Multiple port models equipped with heated wells for cleaning individual bottle with acid vapor

PFA Bottles

- Source and Product bottles are inert to corrosion & reduce contamination

Purification Port(s)

- Some models support multiple vessel kits for simultaneous purification

External Control

- 1.7 m (5.5') PTFE-shielded control cable
- Regulator located remote from CleanAcids module

Graphite Heating Block

- PFA-coated graphite assembly inert to acid corrosion & uniformly heats vessel

Temperature Control

- 110 or 208-240 V supply
- Precise control up to 250°C
- Digital control with programmable timer

Thermal Protection

- PTFE outer shell protects user from heated block

