

# prepFAST Features

## SampleSense Unique Features

### Parabolic Probe

- Perfectly rounded sample transfer line
- Robust, heavy-walled, self-supporting

### Autocorrecting x-y-z

- DXCi autosampler automatically goes to the correct position, every time
- Automatically retries position if an obstacle is encountered

### Highest Speed of Operation

- Corrects for stall events, even after years of use
- Keep analysis running, even after obstruction

## prepFAST Valve Options

### P-Series prepFAST

- High performance sample loading, autodilution and autocalibration

### SampleSense prepFAST

- Automatic optically-sensed sample loading
- Eliminates method timings
- Actively detects and reports sample loading issues

### Dilution Range

- 1 to 400x dilution range (M5)
- 1 to 40,000x dilution range (M5X)

## P-Series and SampleSense Features

### Optional Vented Enclosure

- ULPA filter option

### AutoAlign Arm

- Self-aligns after encountering capped sample tube
- Metal-free fluoropolymer probe aligner

### Sample Containers

- Flexible rack configurations
- Samples from micro vials to large bottles

### Dual Overflowing Rinse Stations

- Ultra-low cross-contamination
- Independently pumped for low rinse consumption

### Vacuum Pump

- Vacuum loaded sampling for rapid uptake

### Chemically Resistant

- Designed for exceptionally long lifetime
- No metallic components above samples

### Syringe Sample Loading

- Precise sample loading for small samples or multiple reanalysis

### Inline Syringe Sample Preparation

- Autocalibration
- Prescribed sample dilution & overage dilution
- Internal standard addition

### Large Volume Reagent Bottles

