

Electronics in Rear Plenum

Optional ULPA Filter

Polypropylene Enclosure

Sample Containers

<1 mL up to 500 mL

Fluoronetic Autosampler

Single Use Resin Columns

Pack and unpack resin automatically

Exhausted Enclosure

Dual Columns

Flow path through single or dual columns

Variable Stir Plates with Resin Bottles

Automatically slurry resin on demand

Syringe Control of Exact Volumes & Flow Rates

For samples and all reagents

Reagent Reservoirs

Exhausted Base Enclosure

