

seaFAST SP2

Electronics in Rear Plenum

Exhausted Enclosure

Built-in to provide maximum contamination control

Sample Containers

<1 mL up to 500 mL

Syringe-Driven Chemistry

Inline or offline, ultrapure preconcentration and matrix removal

Over-Flowing Dual Rinse Station

Minimize sample contamination

Ultraclean Valve System

Completely enclosed liquid flow path

Ultra High Purity Syringe Dilution

Inline, automated calibration and dilution system. Lower blanks in Preconcentration mode.

Bookshelf Reagent Bottles

Maximizes storage for all reagents

