


Clean3V Enclosure & ULPAClean 10

Ultra-clean Sample Environment
Ultra Pure | User-friendly | Versatile



Ultra-clean

- Clean sample environment
- Purifies ambient air (*ULPAClean 10*)
- Protects lab from harmful fumes

User-friendly

- Vertical sliding doors
- Three-sided sample access
- Convenient filter replacement

Versatile

- Variety of sizes
- Available for all ESI autosamplers

4DX Clean3V Enclosure Specifications

	4DX Enclosure	4DX Enclosure with ULPAClean 10
Width	78.7 cm (31 in)	78.7 cm (31 in)
Depth	30.5 cm (12 in)	30.5 cm (12 in)
Height	53.3 cm (21 in)	67.4 cm (26.5 in)
Weight	12.7 kg (28.0 lb)	14.8 kg (32.6 lb)