

New

Heated Rack for ESI's SC-Autosamplers



Heated Rack (5 x 12 Configuration)



Interchangeable Controller

Benefits

- Direct analysis of viscous and semi-solid samples such as base oils, lubricants, vegetable oils, etc. at elevated temperatures
- 2 pre-set temperatures
- Avoids sample prep and diluted solvents
- Controller is transferable between racks
- PTFE insulated
- Optional heated transfer line available
- Optional heated nebulizer and spray chamber

Specifications

- 12V, 120W
- 346.8mm x 105mm x 91.2mm

Temperature

- Standard presets:
 - 0 setting = Ambient
 - 1 setting = 45°C
 - 2 setting = 60°C
- Custom temperature range: 30°C-90°C

Rack Configurations

- 3 x 7 for 30mm o.d. vials (HR-21-30)
- 4 x 10 for 20mm o.d. vials (HR-40-20)
- 5 x 12 for 16mm o.d. vials (HR-60-16)



Heated Rack on SC-2 FAST Autosampler