


MP² Pump Tubing Starter Kit

A wide variety of MP² pump tubing to cover most applications



Kit Includes:

- Full Range of Sizes – 0.13 to 3.17mm i.d. tubing
- Flared and Non-Flared
 - Sizes ≤ 0.76mm are Flared
 - Sizes > 0.76mm are Non-Flared
- PVC, Santoprene and Organic Solvent Resistant
 - 9 PVC Flared tubing packages
 - 17 PVC Non-Flared tubing packages
 - 3 Santoprene Non-Flared tubing packages
 - 4 Solvent Resistant Flared tubing packages
 - 1 Solvent Resistant Non-Flared tubing package



Non-Flared



- 30 CTFE Fluoropolymer Peripump Fittings
 - 5 small male & female fittings for up to 0.89mm i.d. tubing
 - 5 medium male & female fittings for up to 1.75mm i.d. tubing
 - 5 large male & female fittings for up to 3.17mm i.d. tubing

Peripump Fittings			
Part #	ES-2501-PPF2	ES-2501-PPF3	ES-2501-PPF5
Qty	5 Pk	5 Pk	5 Pk
Size	Female (small)	Female (medium)	Female (large)
Part #	ES-2501-PPM2	ES-2501-PPM3	ES-2501-PPM5
Qty	5 Pk	5 Pk	5 Pk
Size	Male (small)	Male (medium)	Male (large)

MPP-K-1 Starter Kit

PVC TUBING

i.d.		Stop Colors		Part #	Qty
0.13mm	Flared	Orange	Black	MPP-013-F-PVC	12 Pk
0.19mm	Flared	Orange	Red	MPP-019-F-PVC	12 Pk
0.27mm	Flared	Orange	Blue	MPP-027-F-PVC	12 Pk
0.38mm	Flared	Orange	Green	MPP-038-F-PVC	12 Pk
0.44mm	Flared	Green	Yellow	MPP-044-F-PVC	12 Pk
0.51mm	Flared	Orange	Yellow	MPP-051-F-PVC	12 Pk
0.57mm	Flared	White	Yellow	MPP-057-F-PVC	12 Pk
0.64mm	Flared	Orange	White	MPP-064-F-PVC	12 Pk
0.76mm	Flared	Black	Black	MPP-076-F-PVC	12 Pk
0.89mm	Non-Flared	Orange	Orange	MPP-089-PVC	12 Pk
0.95mm	Non-Flared	White	Black	MPP-095-PVC	12 Pk
1.02mm	Non-Flared	White	White	MPP-102-PVC	12 Pk
1.09mm	Non-Flared	White	Red	MPP-109-PVC	12 Pk
1.14mm	Non-Flared	Red	Red	MPP-114-PVC	12 Pk
1.22mm	Non-Flared	Red	Grey	MPP-122-PVC	12 Pk
1.30mm	Non-Flared	Grey	Grey	MPP-130-PVC	12 Pk
1.42mm	Non-Flared	Yellow	Yellow	MPP-142-PVC	12 Pk
1.52mm	Non-Flared	Yellow	Blue	MPP-152-PVC	12 Pk
1.65mm	Non-Flared	Blue	Blue	MPP-165-PVC	12 Pk
1.75mm	Non-Flared	Blue	Green	MPP-175-PVC	12 Pk
1.85mm	Non-Flared	Green	Green	MPP-185-PVC	12 Pk
2.06mm	Non-Flared	Purple	Purple	MPP-206-PVC	12 Pk
2.20mm	Non-Flared	Purple	Black	MPP-220-PVC	12 Pk
2.54mm	Non-Flared	Purple	Orange	MPP-254-PVC	12 Pk
2.79mm	Non-Flared	Purple	White	MPP-279-PVC	12 Pk
3.17mm	Non-Flared	Black	White	MPP-317-PVC	12 Pk

SANTOPRENE TUBING (Long life drain tubing)

i.d.		Stop Colors		Part #	Qty
1.30mm	Non-Flared	Grey	Grey	MPP-130-PHR	12 Pk
1.52mm	Non-Flared	Yellow	Blue	MPP-152-PHR	12 Pk
3.17mm	Non-Flared	Black	White	MPP-317-PHR	12 Pk

ORGANIC SOLVENT RESISTANT TUBING

i.d.		Stop Colors		Part #	Qty
0.19mm	Flared	Orange	Red	MPP-019-F-S	12 Pk
0.38mm	Flared	Orange	Green	MPP-038-F-S	12 Pk
0.51mm	Flared	Orange	Yellow	MPP-051-F-S	12 Pk
0.76mm	Flared	Black	Black	MPP-076-F-S	12 Pk
1.30mm	Non-Flared	Grey	Grey	MPP-130-S	12 Pk