

PFA Teflon Evaporation Vessel and Press Fit Cap for Modblock™

Manufactured from chemically inert PFA Teflon, these vessels provide a high purity, low contamination container for performing aggressive acid digestions on the CPI Modblock™ Digestion System. The addition of the o-ring free PFA Purge Port Cap allows the samples to be covered and purged to speed up sample evaporation and reduce contamination from the environment and other nearby samples.

- Directly replaces standard Modblock™ vessels
- PFA-Teflon - High-purity/Low-contamination
- PFA-Teflon resistant to a wide variety of acids
- High temperature (<200°C)
- Samples covered during evaporation
- May be purged with high purity N2 or Ar to speed sample evaporation
- Safely remove acid vapors from Modblock™ area
- Reduce contamination from environment
- Reduce cross-contamination between samples



PFA Evaporation Vessel, PFA Purge Port Cap and PFA storage cap.

Ordering Information	
Part Number	Description
ES-7210	PFA Evaporation Vessel (Pack of 6)
ES-7211	PFA Purge Port Cap (Pack of 6)
ES-7633	PFA Press Fit Cap (Pack of 6)

