

## Precision Quartz Torches and O-ring-free Injectors

Part of a complete o-ring-free sample introduction system

- Direct replacement for standard torches and injectors
- O-ring-free Teflon base
  - reduces contamination
  - easy to clean
  - easy to insert
- Precision alignment of torch and injector eliminates injector wobble and off-axis deflection.
- Precision o-ring-free injector with PFA base available in a variety of sizes and materials including:
  - Quartz
  - Sapphire
  - Platinum
- Stronger, reinforced torch base helps prevent accidental breakage during maintenance
- Available for most ICP and ICPMS instruments



The precision quartz torch and o-ring-free injector directly replace the standard demountable torch and injector.



O-ring-free seal for reduced contamination and easy maintenance.