



prepFAST vs prepFAST S Features Comparison

prepFAST Features

- Autocalibration
- Auto Sample Dilution
- Auto Internal Standard Addition
- Auto QC Overage Dilution



prepFAST

prepFAST S Features

- Autocalibration
- Auto Sample Dilution
- Auto Internal Standard Addition
- Auto QC Overage Dilution
- Automated MSA Calibration
- Automated QC Spiking
- Automated Acid Addition
- Syringe-Driven Sample Introduction
- Enclosed Mobile Stand
- Continuously-Flowing High Purity DI Water Rinse
- DI Water Polishing Columns
- Pre-Cleaned to ppt Level Backgrounds
- PTFE Autosampler Deck



prepFAST S

VS