



Multiple Access Panels

- Left side panel with hinged access
- Optional right side hinged panel
- Front & back panels slide vertically

AutoAlign Arm

- Self-realigns after encountering capped sample tube
- Metal-free fluoropolymer probe aligner

Sample Containers

- <1 mL up to 500 mL

PTFE Autosampler Deck

- Chemically resistant for long lifetime

Position Autocorrection

- Adapts to every position automatically
- Automatically retries position if there is an obstacle
- Keeps analysis running, after years of use

Base Purge Option

- Prevents entry of corrosive fumes into electronics

Optional ULPA Filter or Air Shower

Ultrapure PFA Probe

- Fluoropolymer sample probe

Fluoronetic Drive

- True fluoropolymer mechanical action
 - Magnetically coupled linear drive
- Chemically resistant and long-lived

Ultraclean Enclosure

- Longitudinal PTFE exhaust manifold in deck
 - Fumes removed at sample level

Dual Flowing Ultraclean Rinse Station

- Ultra-low cross-contamination
- Independently pumped for low rinse consumption
 - Flowing UPW rinse option

SampleSense FAST & prepFAST Options

- Automated, ultra-high sample throughput with SampleSense FAST
- Autocalibration, autodilution, and automated sample loading with SampleSense prepFAST