

External Control

- 1.7m (5.5') PTFE-shielded control cable
- Allows regulator to be located outside fume hood

Temperature Regulator

- 110 or 208-240 V supply
- $\pm 2^{\circ}\text{C}$ control up to 250°C maximum
- Multiple alarms for maximum safety
 - Offered in fully programmable and non-programmable (with timer) versions

Two Capacity Options

- 24 vials (29 cm x 19 cm x 21 cm)
- 48 vials (37 cm x 27 cm x 20 cm)

Heating Element

- Embedded in PFA coated graphite block for inertness

Safety

- PTFE thermal insulation around all sides and underneath for safety
- Protective lip surrounds edge of work surface to protect against spills

