

Raised Rim

- Standard raised rim around edge to contain spills

External Control

- 1.7m (5.5') PTFE-shielded control cable
- Allows regulator to be located outside fume hood

Temperature Regulator

- 110 or 208-240 V supply
- $\pm 2^\circ\text{C}$ control up to 250°C maximum
- Multiple alarms for maximum safety
 - Offered in fully programmable and non-programable (with timer) versions

Three Capacity Options

- A4 Plate (small - 20cm x 30cm x 13cm)
- A3 Plate (medium - 30cm x 40cm x 13cm)
- A2 Plate (large - 40cm x 60cm x 13cm)

Heating Element

- Embedded in PFA coated graphite block for inertness

Safety

- PTFE thermal insulation around all sides and underneath for safety
- Protective lip surrounds edge of work surface to protect against spills

