

# prepFAST OilStation Features

## AutoAlign Arm

- Self-aligns after encountering capped sample tube
- Metal-free fluoropolymer probe aligner

## Expanded Sample Deck

- Hundreds of compatible racks and vials
- Compatible with janus racks

## Expanded Storage

- Comfortably fits 20L bottles
- Decrease switching of waste & diluent bottles

## Integrated Touchscreen

- Easy-to-use software
- Stores an unlimited number of custom methods

## 50 mL Syringe Module

- Provides rapid and precise dispensing
- Made from low maintenance, high-quality materials

## Dual Overflowing Rinse Stations

- Ultra-low cross-contamination
- Independently pumped for low rinse consumption

## Integrated Mobile Cart

- Compact system footprint
- Conveniently stores necessary cables

