

# SGTX

## Exhausted Enclosure

- Protects laboratory environment when analyzing concentrated chemicals

## TRAXDisplay

- Easily run methods with LabSymphony software
- Library stores unlimited number of methods

## TRAXFluoroGrip

- Automatically accesses and grips bottles
- Transfers bottles to the TRAXLink for uncapping and barcode reading

## TRAXLink

- Automatic uncapping/recapping
- Barcode reading & sample tracking

## Interlocks

- Automatic UPW shut off
- User alerts for lab safety

## TRAXSpin

- Homogenizes sample prior to analysis
- Mobilizes nanoparticles into solution

## Overflowing UPW Rinse Station

- Rinses probe for baseline blanks


## Instrument Bench

- Metal-free and chemically resistant construction
- Configurable to accommodate most analytical techniques

## Leak Sensors

- External warning light and auto-shut-off

## Consumables Storage



Analytical  
Instrument