

seaFAST SP3



Electronics in Rear Plenum

Exhausted Enclosure

Built-in to provide maximum contamination control

Sample Containers

<1 mL up to 500 mL

Over-Flowing Dual Rinse Station

Minimize sample contamination

Ultraclean Valve System

Completely enclosed liquid flow path

Hydride Generator

Ultratrace detection limits for hydride-forming elements

Syringe-Driven Chemistry

Inline or offline, ultrapure pre-concentration and matrix removal

Ultra High Purity Syringe Dilution

Inline, automated calibration and dilution system. Lower blanks in Preconcentration mode.

Bookshelf Reagent Bottles

Maximizes storage for all reagents

