


ETC EasyTraceCleaner® - Options

Option 1 – Middle Plate

This middle plate doubles your cleaning capacity by allowing you to wash your vessel above your vials.

Middle plate



Doubles your cleaning capacity



Wash your closures above your vials



Designation :	Reference :
Middle Plate for ETC EVO II size M	OPTION-1-M-ETC EVO II
Middle Plate for ETC EVO II size L	OPTION-1-L-ETC EVO II
Middle Plate for ETC EVO II size XL	OPTION-1-XL-ETC EVO II

Option 2 – Additional Tray

This Tray includes 1 lower plate + 1 handle + 48 Rods.


This option allows you to prepare a second set of Vessel to be washed in advance.



Designation :	Reference :
Additional Tray for ETC EVO II size M	OPTION-2-M-120-ETC EVO II
Additional Tray for ETC EVO II size L	OPTION-2-L-120-ETC EVO II
Additional Tray for ETC EVO II size XL	OPTION-2-XL-120-ETC EVO II

Option 3 – 48 Rods


Allows you to stabilize your vessel more easily. Transforms a middle plate for closures into a middle plate for vessel



Designation :	Reference :
Additional 48 Rods in 120mm for ETC EVO II	OPTION-3-120-ETC EVO II
Additional 48 Rods in 160mm for ETC EVO II	OPTION-3-160-ETC EVO II
Additional 48 Rods in 210mm for ETC EVO II	OPTION-3-210-ETC EVO II

Option 4 – P'TIPS RACK

Allows you to easily and efficiently clean all your pipette tips in your ETC.




Designation :	Reference :
P'TIPS Rack 200µl	OPTION-4-200-ETC EVO II
P'TIPS Rack 1000µl	OPTION-4-1000-ETC EVO II
P'TIPS Rack 5000µl	OPTION-4-5000-ETC EVO II
P'TIPS Rack 10000µl	OPTION-4-10000-ETC EVO II
P'TIPS Transfer Rack 200µl	OPTION-4B-TRANSFER-200
P'TIPS Transfer Rack 1000µl	OPTION-4B-TRANSFER-1000

Option 5 – Small OD Vessel Holder

This Vessel holder can be used to stabilize your vessel instead of your rods (picture 2). Essential for holding vessel that cannot be held by the PTFE rods (picture 3).

The height of the holder can be adapted according to the size of the vessel to be cleaned.

Small Vessel Holder



Stabilize your vessel



Hold your vessel without rods



Designation :	Reference :
Small OD Vessel Holder for ETC EVO II	OPTION-5-ETC EVO II

Option 1/5b – Additional Spacer

Allows you to customize the height of your trays (middle plate or Vessel holder) according to your specific needs.

Contact us to define the height of the Spacer in accordance with your Vessel.



Désignation :	Référence :
FOR PTFE MIDDLE PLATE OR VESSEL HOLDER FOR ETC EVO II	OPTION-1/5b-ETC EVO II

Option 6 – Fume Hood for ETC

Helps maintain control of your work environment



Designation :	Reference :
Fume Hood for ETC EVO II	OPTION-7-ETC EVO II

Option 6b – Charcoal Filter

The Charcoal filter system provides optimum safety during your chemical handling



Designation :	Reference :
Charcoal filter for Fume Hood	OPTION-6B-CHARCOAL FILTER