


ULPAclean 10



E-PTFE membrane



ULPAclean 10 is a U15 ULPA filter which is compatible with all DX autosampler enclosures and includes an easily replaceable boron-free, E-PTFE membrane

Ultra Pure

- Purifies ambient air (U15 ULPA filter)
 - Removes 99.9995% of airborne particles with size of 0.1 μm or larger
- E-PTFE filter material
 - Boron Free
- Reduces risk of sample contamination
 - Positively pressurizes enclosure
 - Prevents introduction of airborne contaminants

Convenient Maintenance

- Simple ULPA filter replacement
 - Remove 4 screws, replace filter, replace screws
- Washable pre-filter

Easy to Use

- Two air volume settings
 - 1.0 (m^3/min)
 - 0.5 (m^3/min)
- Quiet operation
- Compatible with all DX autosampler enclosures



Elemental Scientific
ICP | ICPMS | AA