

Physical Details

Dimensions and Weight

- Height: 1544 mm (60.79")
- Width: 794 mm (31.25")
- Depth: 548 mm (21.56")
- Weight: 54.07 kg (119.2 lbs.)

Power Requirements

- Input power: 2 @ 100-240 VAC \pm 10%, 50-60 Hz, 2.8A

Communication

- 2 USB Connections (Instrument computer not included)

Environmental Conditions

- Ambient operating temperature 15°C to 35°C
- Ambient operating humidity 20 to 80% RH non-condensing

Compliance with EMI, Safety Standards

- C.E. Marking: Low Voltage Apparatus 73/23/EEC, EMC89/336/EEC
- Emissions: EN61326, EN55011 Class A, Group 1
- Immunity: EN61326, EN6100-4-2, EN6100-4-3, EN6100-4-4, EN6100-4-5, EN6100-4-6, EN6100-4-11
- The U.S. Code of Federal Regulations, Title 47, FCC Part 15
- Industry Canada ICES-003
- Safety: EN61010-1:2010

8DX prepFAST with ULPA Filtered Enclosure (ULPA Filter and Enclosure are optional)

